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Publisher: Széchenyi István University
Amount of support: 767.08 million Hungarian forints

Produced with the support of the European Agricultural Fund
for Rural Development

Year of publication: 2025.
All rights reserved.

Food Analytical Laboratory

In the laboratory operating at the Albert Kázmér Mosonmagyaróvár Faculty of Széchenyi István University in Győr (SZE-AKMK), researchers conduct small and large-scale instrumental analyses of agricultural and food industry samples.

Sample preparation procedures and small-scale instrumental analyses take place in the balance room, and the inorganic and organic preparation rooms. A separate room is reserved for the large-scale analytical instruments. The equipment and tools for sensory examinations are also located in a separate room.

Check out our videos in the **Smart Garden** app!

Instrumental analytical food examinations

Food safety is one of the most important issues of our time. Every day we encounter news about contaminated, counterfeit, or health-damaging products.

In the food industry, precise, reliable, and rapid analytical testing is becoming increasingly important, as people's health depends on it.

How can the nutritional value of a food be determined? How do we know how many calories we consume?

Analytical tests help uncover the true composition, quality, and potential contaminants of foods.

Main examinations in the analytical laboratory

- ◆ dry matter and ash content
- ◆ determination of crude protein, fat content, and fatty acid composition
- ◆ pH measurement
- ◆ determination of total polyphenol and total antioxidant content
- ◆ examination of organic acid and sugar composition
- ◆ the effect of different drying methods (lyophilisation, drying, oven drying) on the nutritional value of foods

The analytical laboratory at Széchenyi István University is at the forefront of applying the most modern techniques and plays a significant role in research and development.

Gravimetric and titrimetric methods

These classical methods of chemical analysis provide relatively quick and reliable results during examinations. The gravimetric method is based on mass measurement, while the titrimetric method is based on volume measurement.

During the **gravimetric method**, the mass of a specific substance is measured after drying, ashing, destruction, or solvent extraction, and the result is compared to the mass of the original (untreated) sample. For example, this method can determine how much ash remains after burning a food – this can indicate the mineral content. The essence of **titrimetry** is that a solution of known concentration (titrant) is dripped into the sample being tested until a specific, easily detectable chemical reaction (change in colour, pH, or electromotive force) occurs. In this way, for example, the acidity of a fruit juice can be precisely determined.

Applications of titrimetry for determining:

- ◆ acidity (e.g. vinegar, wine),
- ◆ amount of vitamin C,
- ◆ the amount of calcium in milk.

Smart Garden app: the method of „home titration”.

Spectrophotometric methods

Spectrophotometry is a much more advanced procedure. The principle is simple: different substances (atoms, molecules) absorb or emit light at different wavelengths during their energy changes. If a substance is illuminated with light of a given wavelength, its absorption can tell us what chemical components are present in the sample and in what quantities.

This method is suitable for determining the colouring agent content of foods, estimating protein content, and detecting nitrate and nitrite content.


Smart Garden app: detecting honey adulteration using the refractometry method.

Chromatographic methods

Chromatographic methods are among the most advanced techniques. Chromatography is a physical separation method used to separate and identify components of complex samples. Its essence is that the components to be separated are distributed between two phases, a stationary phase, and a mobile phase moving in a definite direction. The method is based on physicochemical properties (e.g. solubility, adsorption, mass).

General structure of a gas chromatograph

The **HPLC systems** are capable of detecting even **nanogram** quantities – it's like searching for a grain of sugar in a swimming pool full of water!


Chromatography – main types

- ♦ **Gas chromatography (GC):** for examining volatile substances, aromas, contaminants

Gas chromatography is a forced-flow separation technique where the mobile phase is in a gaseous state. Since we use a gaseous mobile phase to flow the components, they must also be in a gaseous state. It follows that gas chromatography is fundamentally used for the analysis of volatile components.

- ♦ **High-performance liquid chromatography (HPLC):** for detecting non-volatile, heat-sensitive substances (such as vitamins, sweeteners)

Modern HPLC device and the heart of the device (stationary phase or column – where the separation occurs)

Development of the laboratory at Széchenyi István University

The profile of the Food Analytical Laboratory includes a variety of methods from simple gravimetric and titrimetric tests, through spectrophotometric tests, to large-scale chromatographic tests.

The main research areas of the laboratory:

- ◆ Development and optimisation of extraction methods,
- ◆ Examination of the total polyphenol content and antioxidant effect of plant extracts, fruits, and vegetables,
- ◆ Development and spectrophotometric examination of foods supplemented with plants,
- ◆ Oil extraction from seeds using Soxhlet extraction,
- ◆ Examination of sugars and organic acids, phenolic acids,
- ◆ Extraction and composition analysis of plant essential oils,
- ◆ Examination of the sugar, organic acid, fibre, and volatile acid content of fermented feeds,
- ◆ Examination of sweeteners in soft drinks.

In addition to research projects, the laboratory plays a significant role in the practical education of students in food engineering (BSc) and food safety and quality engineering (MSc).

Within the framework of its grant-winning rural development project, the University aims to further develop the institution's cheese and fruit processing plant production using herb and spice crops grown in the newly established herb garden on the Mosonmagyaróvár campus.

The industrial applicability of the herbs intended for use in the small-scale dairy processing plant can be monitored partly by changes in microbiological and compositional parameters (such as acidification properties) and partly by organoleptic properties. Determining the compositional parameters of raw materials, as well as the compositional parameters of semi-finished and finished products, is nowadays primarily based on Fourier transform infrared technology, applied in a non-contact manner for finished products, alongside standard, conventional methods. To demonstrate the technology, two devices were procured in the project.

The two new devices are necessary because the medium to be examined is a liquid sample matrix in the case of raw milk, while it is a solid sample matrix for semi-finished and finished products. In the former case, the device assisting the examinations is the **MilkoScan**, and in the latter, the **FoodScan**.

Among the raw materials used for the herbal products developed within the project, determining the essential oil content and, within that, the **active ingredient content** of the **herb and spice plants** is of paramount importance. This is demonstrated by the new device procured within the project, based on an analytical method: the gas chromatograph-mass spectrometer (GC-MS).

Determination of essential oil and active ingredient content using GC-MS

The active ingredient content and physiological properties of herbs and spices largely depend on the composition of their essential oils, as these contain the aroma compounds and biologically active components. The GC-MS procedure allows for their identification and quantitative determination. Steps of the examination:

Sample preparation – extraction of essential oil

The aim of sample preparation is to separate and concentrate the compounds to be determined from the sample matrix composed of unwanted components.

Smart Garden app: video about the GC-MS device.

Gas Chromatographic Separation (GC)

The sample is injected into the GC device's injector, where due to the high temperature (200–300 °C), the sample components become gaseous, then proceed through a long column.

Since individual molecules travel at different speeds in the column due to their varying chemical properties (e.g., polarity), they **separate from each other**.

Mass Spectrometric Identification (MS)

The molecules exiting the column are ionised by the mass spectrometer and analysed according to their **mass/charge ratio**. This allows for precise determination of what compounds are present, even if they are only in very small quantities.

Evaluation – creating an essential oil profile

The result of the GC-MS measurement is a detailed essential oil map. Based on the obtained chromatogram (peak diagram), it can be determined: which essential oils are present (identification based on a database), and in what proportion (based on peak area).

With the help of the device, it can be easily determined whether a given product (e.g., lavender oil) truly originates from a herb or is just a scented imitation.

Dairy raw material and product testing with MilkoScan and FoodScan devices

During the production of dairy products, it is important to quickly obtain information about the composition of the raw material and the product to ensure compliance with regulatory criteria, optimise quality, and maximise product yield, thereby increasing profit.

In **butter production**, it is a legal requirement that the product is made with a minimum of 80% fat content or a maximum of 16% moisture content.

The most effective rapid testing methods for determining the compositional values of industrial and laboratory dairy products are near-infrared (NIR) spectroscopy and Fourier-transform infrared (FTIR) spectroscopy.

The **MilkoScan™** Mars is an easy-to-use and cost-effective

device, which, through **Fourier-transform infrared technology**, enables the laboratory to measure six parameters: fat, protein, lactose, total solids, and

fat-free dry matter content, and exclusively for milk – to detect adulteration – freezing point, from raw and processed milk, as well as whey and cream, without the


need for complex sample preparation procedures, with a measurement time of 1 minute.

FoodScan™ 2 Dairy enables the rapid analysis of dairy products and plant-based dairy imitations through **near-infrared**

technology. Instead of a single examination, the device scans the sample three times, thus examining a larger portion of the sample, making the result more representative of the compositional values characteristic of the entire sample mass.

What does the FoodScan measure?

- ◆ fat and moisture content
- ◆ protein and salt content
- ◆ fat-free dry matter
- ◆ total solids content
- ◆ the saturated fatty acid content of cheeses

You can read about sensory food testing in our **Smart Garden** app.

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